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rProtein A HC-650F 亲和层析填料

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HIGH
QUALITY
EXPERIENCE

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
3. 使用 0.1 N NaOH 对 TOYOPEARL AF-rProtein A HC-650F 在位清洗 (CIP)


层析柱: TOYOPEARL AF-rProtein A HC-650F (5 mm I.D. x 5 cm)

CIP 循环: 10 CV with binding buffer, 5 CV with elution buffer, 7 CV with binding buffer

3 CV with 0.1 N NaOH (contact time; 15 min), 5 CV with binding buffer

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

样品：40 mL feedstock from cell culture (humanized IgG1: 1 g/L)

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Product name	Supplier	Base matrix	Particle size (micron)	Binding capacity (g/L)*			Chemical stability	List price (\$/L)**	Comment
				SBC	DBC (2 min)	DBC (5 min)			
TOYOPEARL AF-rProtein A HC-650M	Tosoh	Polymer	45	>68	50	70	NaOH durable	NA	rProtein A; C-domains, Multipoint attachment
TOYOPEARL AF-rProtein A 650M	Tosoh	Polymer	45	>45	30	40	NaOH durable	12,000	rProtein A; 4 C-domains, Multipoint attachment, high flow type
MabSelect SuRe LX	GE	Agarose	85	ND	30	58	NaOH durable	16,000	rProtein A; 4 B-domains, Single point attachment
MabSelect SuRe	GE	Agarose	85	ND	30	40	NaOH durable	15,000	rProtein A; 4 B-domains, Single point attachment
AbSolute	NovaSep	Silica glass	44	70	40	45	pH 1.5 - 9.0	18,000	High flow type, Bed height max 40 cm, 1 % benzyl alcohol/acetate buffer, < 10 CFU/mL
ProSep UltraPlus	Merck	Porous Glass	60	> 67	45	50	pH 1.5 - 8.5	15,000	rProtein A; natural type, Multipoint attachment, High flow type, Bed height max 30 cm, Shipped with 1 % benzyl alcohol/acetate buffer
KanCapPro A	Kaneka	Cellulose	75	ND	30	50	NaOH durable	?	rProtein A; 5 C-domains, announced to expand production scale in 2013
Amsphere Protein A JWT203	JSR	Polymer	50	ND	33	43	NaOH durable	?	Hydrophobic properties?

* Catalogue data
** Estimated price

TOSOH CORPORATION, BIOSIENCE DIVISION

- GE MabSelect SuRe LX
- Novasep AbSolute
- Merk ProSep UltraPlus
- Thermo/Life Technologies, POROS Protein A
- Kaneka KanCap Pro A

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* ELISA Kit is different from that for TOYOPEARL AF-rProtein A 650F.

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